

CellCounter

Fetal Hgb Counter

Cell counter for measuring the fetal hemoglobin percentage in a maternal sample. Also known as the Kleihauer-Betke (KB) test

Step 1:

Count the number of RBCs in 5 fields using the high dry objective (40X);

Total RBCs in 5 fields:	Average RBCs per field:	RBCs in 30 fields:
21	4.2	126

Field 1:	Field 2:	Field 3:	Field 4:
Total RBCs (Adult and Fetal) 6	Total RBCs (Adult and Fetal) 4	Total RBCs (Adult and Fetal) 4	Total RBCs (Adult and Fetal) 4

Field 5:
Total RBCs (Adult and Fetal) 3

Step 2:

Count the number of fetal cells in 30 fields using the high dry objective (40X);

Percent Fetal cells:

3.2

Fetal Cells in 30 Fields:	
Fetal Cells: 4	Fields: 3

Notice: CellCounter is intended to be used as a visual tallying aid for manual cell counting. Any calculations displayed by the tool are provided

for informational and confirmatory purposes only and must be independently verified in accordance with institutional SOPs prior to result reporting.

Wednesday, January 7, 2026 at 21:25:43

© 2024-2026 CellCounter. All rights reserved. | Version 1.0

[Terms of Use](#)