

CellCounter

Hemocytometer Counter

Cell counter for tallying the amount of red blood cells and white blood cells in a body fluid.

Final RBC Count:

0

Final TNC Count:

0

Squares Counted: 9

Dilution Factor: 1

Side 1

RBC

0

TNC

0

Side 2

RBC

0

TNC

0

☐ Click here to use separate dilutions and square counts for RBC and TNC

☐ Click here to show preliminary calculations

RBC Count = $\frac{(\text{RBC Average}) \times (\text{dilution}) \times 10}{(\text{\# of Squares Counted})} = \frac{(0) \times (1) \times 10}{(9)} \approx 0$

TNC Count = $\frac{(\text{TNC Average}) \times (\text{dilution}) \times 10}{(\text{\# of Squares Counted})} = \frac{(0) \times (1) \times 10}{(9)} \approx 0$

Notice: CellCounter is intended to be used as a visual tallying aid for manual cell counting. Any calculations displayed by the tool are provided for informational and confirmatory purposes only and must be independently verified in accordance with institutional SOPs prior to result reporting.

Wednesday, January 7, 2026 at 21:20:42

© 2024-2026 CellCounter. All rights reserved. | Version 1.0

[Terms of Use](#)